


Circuito odontológico Actynox® Aspiração Passiva

Conexão em espiga na garrafa



4

3

2

8

8

9

5

6

7

10

1A

Máscara de Adultos

1B

Máscara Pediátrica

4

11














12

Tubo corrugado
Ø22 saída
exterior
(NÃO INCLUÍDO)














Montar no sentido do gás conforme indicado



6014460 - Circuito odontológico Actynox® Aspiração Passiva - Adultos - Componentes - Lote: xxxxxx

| Item | Cod. | Lot. | Descrição |  |  | Quant. Circuito |
|------|-----------|------|---|---|---|-----------------|
| 1A | SIL-MED-1 | |  <p>Máscara Odontológica Silhouette</p> <p>CE 0413</p> | | <p>Parker Hannifin Corp. 245 Township Line Road, Hatfield, Pennsylvania 19440, United States</p> | 5 |
| 2 | 2612001 | |  <p>Dispositivo standard (Inclui uma bolsa reservatório de 2l e tubo de condução de oxigénio de 1,8m)</p> <p>CE 0120</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 1 |
| 3 | 1971000 | |  <p>Conetor 22M - 30F</p> <p>CE 0120</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 1 |
| 4 | 1568000 | |  <p>Conetor reto 22F - 6MM</p> <p>CE 0120</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 2 |
| 5 | 4222/701 | |  <p>Filtro bacteriano / vírico</p> <p>CE 2797</p> | | <p>GVS Filter Technology UK Limited NFC House Vickers Industrial Estate Mellishaw Lane Morecambe Lancashire LA3 3EN United Kingdom</p> | 5 |
| 6 | M143 | |  <p>Conetor de Máscara Silhouette da linha inspiratória</p> <p>CE</p> | | <p>Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España</p> | 1 |
| 7 | M144 | |  <p>Conetor de Máscara Silhouette da linha expiratória</p> <p>CE</p> | | <p>Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España</p> | 1 |
| 8 | MV019 | |  <p>Linha de condução de oxigénio 1,8 m</p> <p>CE 0123</p> | | <p>SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China</p> | 2 |
| 9 | 00100510 | |  <p>Conexão "Y"</p> <p>CE 0843</p> | | <p>QUIROSA, S.A. Crta. Artés a Manresa, km 1,5 P.I. Sta. Mª Artés, Nau 16-17 08271 ARTÉS (Barcelona) Spain</p> | 1 |
| 10 | MV003 | |  <p>Linha de condução de oxigénio 10 cm</p> <p>CE 0123</p> | | <p>SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China</p> | 1 |
| 11 | 1950000 | |  <p>Válvula unidirecional</p> <p>CE 0120</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 1 |

6014461 - Circuito odontológico Actynox® Aspiración Passiva – Pediátrico - Componentes - Lote: xxxxxx

| Item | Cod. | Lot. | Descrição |  |  | Quant. Circuito |
|------|----------|------|--|---|---|-----------------|
| 1B | SIL-SM-1 | |  <p>CE 0413</p> <p>Máscara Odontológica Silhouette</p> | | <p>Parker Hannifin Corp. 245 Township Line Road, Hatfield, Pennsylvania 19440, United States</p> | 5 |
| 2 | 2612001 | |  <p>CE 0120</p> <p>Dispositivo standard (Inclui bolsa de reservatório de 2l e tubo de condução de oxigénio de 1,8m)</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 1 |
| 3 | 1971000 | |  <p>CE 0120</p> <p>Conetor 22M - 30F</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 1 |
| 4 | 1568000 | |  <p>CE 0120</p> <p>Conetor reto 22F - 6MM</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 2 |
| 5 | 4222/701 | |  <p>CE 2797</p> <p>Filtro bacteriano / vírico</p> | | <p>GVS Filter Technology UK Limited NFC House Vickers Industrial Estate Mellishaw Lane Morecambe Lancashire LA3 3EN United Kingdom</p> | 5 |
| 6 | M143 | |  <p>CE</p> <p>Conetor de Máscara Silhouette da linha inspiratória</p> | | <p>Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España</p> | 1 |
| 7 | M144 | |  <p>CE</p> <p>Conetor de Máscara Silhouette da linha expiratória</p> | | <p>Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España</p> | 1 |
| 8 | MV019 | |  <p>CE 0123</p> <p>Linha de condução de oxigénio 1,8 m</p> | | <p>SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China</p> | 2 |
| 9 | 00100510 | |  <p>CE 0843</p> <p>Conexão "Y"</p> | | <p>QUIROSA, S.A. Crta. Artés a Manresa, km 1,5 P.I. Sta. Mª Artés, Nau 16-17 08271 ARTÉS (Barcelona) Spain</p> | 1 |
| 10 | MV003 | |  <p>CE 0123</p> <p>Linha de condução de oxigénio 10 cm</p> | | <p>SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China</p> | 1 |
| 11 | 1950000 | |  <p>CE 0120</p> <p>Válvula unidirecional</p> | | <p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p> | 1 |